

FLOWFAST H

TECHNICAL DATASHEET



TECHNICAL SPECIFICATIONS

PRODUCT PROTECTION WORKBENCH

FlowFAST H | Class III according to DIN EN ISO 14644-1

Description	Unit	0.9 m	1.2 m	1.5 m	1.8 m
General dimensions ¹ W×H×D	mm	945 × 1195 × 925	1250 × 1195 × 925	1555 × 1195 × 925	1860 × 1195 × 925
Usable dimensions ² W×H×D	mm	885 × 734 × 600	1190 × 734 × 600	1495 × 734 × 600	1800 × 734 × 600
Unit height with underframe (standard)	mm	2025			
Working height (seated workstation)	mm	880			
Work opening	mm	735			
Weight	kg	104	122	171	194
Laminar flow	m/s	0,40			
Noise level	dB(A)	49	51	53	55
Light intensity	lux	>1100	>1200	>1300	>1300
Electrical connection	1Ph+E - 230V 50Hz				
Power consumption	W	320	460	460	520
Protection class / IP	1/20				
Internal protection	A	4			
Sockets IP66		1	1	1	1
Additional sockets		●	●	●	●
PVC front window		●	●	●	●
Automatic flow regulator		●	●	●	●

¹ The devices are delivered assembled, but can be disassembled on site for installation. Reduction of the total depth to 700 mm possible.

² The devices are also available with a usable device depth of 1000 and 1500 mm.

● Standard ● Optional