



FLOWFAST H REVERSE

TECHNICAL DATASHEET




MADE IN ITALY

TECHNICAL SPECIFICATIONS

PRODUCT PROTECTION WORKBENCH

FlowFAST H Reverse | Class III according to EN ISO 14644-1

Description	Unit	0.9 m	1.2 m	1.5 m	1.8 m
General dimensions ¹ W×H×D	mm	945 × 1195 × 925	1250 × 1195 × 925	1555 × 1195 × 925	1860 × 1195 × 925
Usable dimensions ² W×H×D	mm	885 × 734 × 600	1190 × 734 × 600	1495 × 734 × 600	1800 × 734 × 600
Unit height with underframe (standard)	mm	2025			
Working height (seated workstation)	mm	880			
Work opening	mm	735			
Complete front window opening	mm	480			
Weight	kg	104	122	171	194
Inflow	m/s	>0,5			
Noise level	dB(A)	49	51	53	55
Light intensity	lux	>1100	>1200	>1300	>1300
Electrical connection	1Ph+E - 230V 50Hz				
Power consumption	W	320	460	460	520
Protection class / IP	1/20				
Internal protection	A	4			
Sockets IP66		1	1	1	1
Additional sockets		●	●	●	●
PVC front window		●	●	●	●
Automatic flow regulator		●	●	●	●

¹ The devices are delivered assembled, but can be disassembled on site for installation. Reduction of the total depth to 700 mm possible.

² The devices are also available with a usable device depth of 1000 and 1500 mm.

● Standard ● Optional