



FLOWFAST V

TECHNICAL DATASHEET




MADE IN ITALY

TECHNICAL SPECIFICATIONS

PRODUCT PROTECTION WORKBENCH

FlowFAST V | Class III according to DIN EN ISO 14644-1

Description	Unit	0.9 m	1.2 m	1.5 m	1.8 m
General dimensions W×H×D	mm	1045 × 1545 × 760	1350 × 1545 × 760	1655 × 1545 × 760	1960 × 1545 × 760
Usable dimensions W×H×D	mm	893 × 755 × 609	1198 × 755 × 609	1503 × 755 × 609	1808 × 755 × 609
Unit height with underframe (standard)	mm	2345			
Working height (seated workstation)	mm	865			
Work opening	mm	250			
Complete front window opening	mm	480			
Weight	kg	160	192	274	325
Downflow	m/s	0,45			
Noise level	dB(A)	<52	<55	<55	<55
Light intensity	lux	>1200	>1200	>1300	>1300
Electrical connection	1Ph+E - 230V 50Hz				
Power consumption	W	300	450	550	650
Protection class / IP	1/20				
Internal protection	A	4			
Sockets IP66		1	1	2	2
Additional sockets		●	●	●	●
UV lamp including aluminum cover		●	●	●	●
Automatic flow regulator		●	●	●	●

● Standard ● Optional